

## **Blood Required for Adult Tubes for Different Tests and Test Combinations**

# Category A (1/3 tube, or ∼1.5 mL per bullet point)



### Single analytes

- C-Reactive Protein (CRP)
- Beta-hCG (βhCG)
- Estradiol (E2)
- Insulin
- Luteinizing Hormone (LH)
- Procalcitonin



### **Different test combinations**

- Repro Endo Profile 1 +
   1 Category B test
- Repro Endo Profile 2 +
   1 Category B test
- HIV Screen + HBsAg + Anti-HBs
   + HCV Total Antibody + Syphilis

# Category B (1/2 tube, or ~2.5 mL per bullet point)



### Multiple analytes (Profile test / Panel test)

- Renal Panel (RP)
- Liver Function Test (LFT)
- Lipid Panel
- Iron Panel
- Thyroid Panel (FT4 + TSH)
- Hepatitis B Immunity Profile
- Repro Endo Profile 1
- Repro Endo Profile 2
- Free Testo (Calculation) + FAI
- Cardiac Markers (CK + CKMB + Troponin I)
- Antenatal (HBsAg + Syphilis + HIV Screen + Rubella IgG)

- AFP + βhCG
- FSH + LH
- FSH + LH + E2
- Prolactin + PEG
- Vitamin B12 + Folate +
   Iron Panel
- HBsAg + Syphilis + Anti-HCV
  - + HIV Screen

- CRP + Procalcitonin
- RP + CRP + Procalcitonin
- LFT + RP + Calcium + CEA
- LFT + RP + CRP + Calcium + Magnesium + Phosphate + Uric Acid

### **DISCLAIMER:**

- For tests not listed in this guide, please refer to eLab book.
- Sharing of tubes across different laboratories is not recommended. Please send separate tubes accordingly.
- Blood volume is calculated based on 50% PCV / HCT, and minimum serum height required according to manufacturer's criteria.
   Note: volume of blood required for pediatric tube is different due to smaller tube diameter.
- If the patient has higher PCV / HCT (>50%), it may not be sufficient for half tube to run Category B tests. It is recommended to draw full tube of blood.
- For samples with low volume (i.e. patients who are difficult to bleed), please send a memo to lab. Failure to do so will result in sample rejection.

  Refer to eLab book > Special Instructions > Clinical Chemistry